



Student Iris Scissors - Curved/Sharp-Sharp/Student Stainless Steel/11.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91461-11
Category	Student Line > Scissors & Micro Scissors Most Popular > Student Line
EAN	4251678878156



Find 91461-11 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		