



Student Iris Scissors - Curved/Sharp-Sharp/Student Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

This product is discontinued and has limited availability.

General Info

Product Number	91461-12
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510602



Find 91461-12 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	12 cm		
Feature	Student Quality		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Curved		