



Student Wagner Scissors - Straight/Sharp-Sharp/Student Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91468-12
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510589



Find 91468-12 on finescience.com

Technical Info

Shipping Weight	19 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	30 mm		
Cuttingedge	20 mm		
Length	12 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		