



Student Spencer Ligature Scissors - Straight/Blunt-Blunt/Student
Stainless Steel/9cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91476-09
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510619



Find 91476-09 on
finescience.com

Technical Info

Shipping Weight	14 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Feature	Student Quality		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		