



Student Vannas Spring Scissors - Straight/Sharp/Student Stainless Steel/9cm/5mm Cutting Edge



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91500-09
Category	Student Line > Scissors & Micro Scissors Most Popular > Student Line
EAN	4251909508791



Find 91500-09 on finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	9 cm		
Feature	Student Quality		
Serrated	No		
Tipdiameter	0.35 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		