



Student Spring Loaded Retractor - Blades/Student Stainless Steel/8.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91706-05
Category	Student Line > Retractors
EAN	4251909508845



Find 91706-05 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	5 cm		
Feature	Student Quality		
Maximumspread	3 cm		
Teeth	Atraumatic blades		
Teethdepth	14 mm		
Teethwidth	1 mm		