



Student Senn-Miller Retractor - 3 blunt teeth 7 x 8 mm and 1 blunt tooth 18 x 6 mm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91791-16
Category	Student Line > Retractors
EAN	4251909511241



Find 91791-16 on finescience.com

Technical Info

Shipping Weight	27.5 g	Net Weight	26.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	16 cm		
Feature	Student Quality		
Teeth	3 blunt teeth 7 x 8 mm and 1 blunt tooth 18 x 6 mm		
Teethdepth	5.5 mm		
Teethwidth	2.5 mm		
Tipdimensions	10.5 x 8.5 mm		