



Student Dieffenbach Clamp - Straight/35mm

Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.



General Info

Product Number	91850-35
Category	Student Line > Clamps & Hemostats
EAN	4251909510855



Find 91850-35 on
finescience.com

Technical Info

Shipping Weight	4 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clamppress	Strong		
Length	35 mm		
Feature	Student Quality		
Tipdimensions	13 x 4 mm		
Tipwidth	4 mm		
Tips	Serrated		
Tipshape	Straight		