



Ring Forceps - with Screw/Smooth/Straight/9cm/ID: 4.8 mm - OD: 6 mm



For transport of tissue and organs. A limit screw prevents excessive applied pressure.

General Info

Product Number	11106-09
Category	Forceps & Tweezers > Smooth
EAN	4251909501938



Find 11106-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	9 cm		
Tipdimensions	ID: 4.8 mm OD: 6 mm x 1.5 mm		
Tips	Ring		
Tipshape	Straight		