



Dumont #5AC Forceps - Biology Tips/Straight/Inox/11.5cm



Allgemein

Produktnummer	11258-20
Kategorie	Forceps & Tweezers > Smooth
EAN	4251909502287



Finde 11258-20 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	9 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	Type #5 Anti-Capillary		
Tipdimensions	0.07 x 0.02 mm		
Tips	Biology		
Tipshape	Type #5 Anti-Capillary		