



Kelly Hemostat - Serrated/Curved/14cm



Allgemein

Produktnummer	13019-14
Kategorie	Clamps & Clips > Hemostats
EAN	4251909503369



Finde 13019-14 auf
finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	25 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	22 mm		
Length	14 cm		
Tipwidth	1.5 mm		
Tips	Serrated		
Tipshape	Curved		