



Wired Liebermann Retractor - 7cm/4cm Maximum Spread



General Info

Product Number	17010-07
Category	Retractors
EAN	4251909505929



Find 17010-07 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	8 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	7 cm		
Feature	Wire		
Maximumspread	4 cm		
Shape	Slightly curved down		
Teeth	None		
Teethdepth	3 mm		
Teethlength	4 mm		