



Biemer Clip Applicator Forceps Style - 14cm/Straight/without Lock



Application forceps for Biemer clamps made of surgical stainless steel. Required to use our items FST # 18058-15 & 18058-18.

General Info

Product Number	18059-14
Category	Clamps & Clips > Appliers MedTech Line
EAN	4251909511692



Find 18059-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	38 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	14 cm		
Lockable	No		
Sterile	No		