Pre-treatment Gel Spray



NEW ITEM

Use the pre-treatment gel spray to maintain moisture and protect instruments during transport. This spray can be used at Point-of-Use when delays are likely to occur and when instruments are transferred to instrument reprocessing. Use in final Central Sterile instrument processing area for spot treatment of heavy residual soils or when high case-loads lead to delays. Use on the most challenging orthopedic, neuro-spinal, gastrointestinal cases. The solution is a ready-to-use neutral pH gel spray designed for pre-treating surgical instruments and medical devices. Humectant gel helps prevent soils from drying and hardening after instrument use. Maintains moisture on contaminated instruments for up to 72 hours. Protects delicate surfaces from corrosion during prolonged hold times. Improves cleaning effectiveness and efficiency during subsequent cleaning steps. Formulated with biodegradable surfactants and wetting agents. Compatible with a wide range of materials, including stainless steel, aluminum, plastics, glass and delicate ophthalmic instruments. Benefits include no sticky residue left on carts or trays, rinses off easily during cleaning, replaces the need for messy liquid soak pans, bag-on-valve technology provides broader spray coverage with controlled delivery, and the green-tinted gel foam allows visible confirmation of coverage for compliance.

General Info

Product Number	29040-01
Category	New
EAN	0

Find 29040-01 on
finescience.com

Technical Info

Shipping Weight	612 g	**Net Weight**	500 g
Package Size	1 Pc.	**Unit of Measure**	EA
Minimum Order Quantity	1		
Volume	355 ml		